



Diagnostic Instruments



CROWN CLICK

INTERNATIONAL

CERTIFICATE



Dynamex Certification Services certify that
CROWN CLICK INTERNATIONAL

Registration No. DCS/2232-A

Address to which this certificate refers:

Mehran Street, Mohallah Prem Nagar, Sialkot, (51310) – Pakistan.

Approved Scope to which this certificate refers:

“Manufacturer and Exporter of Surgical & Dental Instruments”

has been evaluated and found to meet the requirements of

ISO 9001: 2015

CHAIRMAN

SCHEME MANAGER

Certificate Issue Date: 08 August, 2024


Certificate Expiry Date: 07 August, 2025

Date of Initial Registration: 08 August, 2024

Re-assessment Date: 07 August, 2027

This Certificate of Registration is granted subject to the Regulations approved by the Board

www.dynamexcertification.org



CE Compliance Certificate



Application of Medical Device Regulation MDR 2017/ 745 for Class I Medical Devices.

This is to certify that the products submitted are:

**MEDICAL DEVICES CLASS I
(Re-Useable Surgical and Dental Instruments)
Registration no DCS/1248322558-A**

Manufactured By:

CROWN CLICK INTERNATIONAL

Mehran Street, Mohallah Prem Nagar, Sialkot, 51310 – Pakistan.

Comply with the applicable requirements of the MDR 2017/ 745, The Technical file of the devices have been assessed, including the EC declaration of conformity according to the Annex IV confirming that the medical devices are fulfilling the applicable requirements of MDR 2017/745.

Limitations:

The manufacturer must inform DCS of any substantial changes occurred in the Product or process in order to examine whether this certificate remains valid.

CHAIRMAN


SCHEME MANAGER

Certificate Issue Date: August 08, 2024

Certificate Expiry Date: August 07, 2025

This Certificate of Registration is granted subject to the Regulations approved by the Board

www.dynamexcertification.org



CERTIFICATE



Dynamex Certification Services certify that
CROWN CLICK INTERNATIONAL

Registration No. DCS/2232-B

Address to which this certificate refers:

Mehran Street, Mohallah Prem Nagar, Sialkot, (51310) – Pakistan.

Approved Scope to which this certificate refers:

“Manufacturer and Exporter of Surgical & Dental Instruments”

has been evaluated and found to meet the requirements of

ISO 13485: 2016

CHAIRMAN

SCHEME MANAGER

Certificate Issue Date: 08 August, 2024


Certificate Expiry Date: 07 August, 2025

Date of Initial Registration: 08 August, 2024

Re-assessment Date: 07 August, 2027

This Certificate of Registration is granted subject to the Regulations approved by the Board

www.dynamexcertification.org





CROWN CLICK INTERNATIONAL



"Businesses require 'Continuous Professional Development' at each activity process to ensure quality growth in its operations. We utilise the years of experience with assistance of the information technology to evaluate and improve year by year. Our aim is to supply our customers with best quality and to bring revolutionary solutions to healthcare industry. We are confident that with our carefully selected team, and their hardwork, we would be able to reach our target and become the best possible source for healthcare equipment from Sialkot, Pakistan."

*Muhammad Ahsan CEO
Crown Click International*

Diagnostic Instruments

for ENT & Ophthalmic Examination

Diagnostic Instruments Set

Full featured 2.5V set with 5 different instruments packed in hard plastic case

02-3050

Ophthalmoscope Head

- is equipped with 29 corrective lenses from -35 to 40 diopter

Otoscope Head

- Has fitting for insufflator
- Comes with removable 3x magnifying lens
- Ear speculum/cannula with internal diameter \varnothing 2, 3 & 4 mm, reusable, made from polycarbonate, autoclaveable at 121° C

Battery Handle

- Ensure connection to instrument heads with a bayonet locking system
- Require 2 "C" type batteries, not included
- Has knurled finish for better grip
- Has a rheostat to vary light intensity
- Spare lamp

Nasal Speculum

- Connects to the otoscope head and comes with manual screw for opening and closing of speculum

Bent Arm Illuminator

- Light the oral cavity and pharyngeal area
- Attaches to laryngeal mirrors No. 2 & 4 (included)
- Attaches to one metal tongue depressor (included)
- Attaches to one reusable adaptor (included) that holds disposable tongue spatula



Otoscope

Packed in hard plastic case

02-3010

Otoscope Head

- Has fitting for insufflator
- Comes with removable 3x magnifying lens
- Ear speculum/cannula with internal diameter \varnothing 2, 3 & 4 mm, reusable, made from polycarbonate, autoclaveable at 121° C

Battery Handle

- Ensures connection to instrument heads with a bayonet locking system
- Requires 2 "C" type batteries, (not included)
- Has knurled finish for better grip
- Has a rheostat to vary light intensity
- Spare lamp



Otoscope / Ophthalmoscope Combination Set

Packed in hard plastic case

02-3030

Otoscope Head

- Has fitting for insufflator
- Comes with removable 3x magnifying lens
- Ear speculum/cannula with internal diameter Ø 2, 3 & 4 mm, reusable, made from polycarbonate, autoclaveable at 121° C

Ophthalmoscope Head

- is equipped with 29 corrective lenses from -35 to 40 diopter

Battery Handle

- Ensures connection to instrument heads with a bayonet locking system
- Requires 2 "C" type batteries, not included
- Has knurled finish for better grip
- Has a rheostat to vary light intensity
- Spare lamp



Ophthalmoscope

Packed in hard plastic case

02-3020

Ophthalmoscope Head

- is equipped with 29 corrective lenses from -35 to 40 diopter

Battery Handle

- Ensures connection to instrument heads with a bayonet locking system
- Requires 2 "C" type batteries, not included
- Has knurled finish for better grip
- Has a rheostat to vary light intensity
- Spare lamp



Diagnostic Instruments

for ENT & Ophthalmic Examination

Mini Otoscope / Ophthalmoscope Combination Set Packed in hard plastic case

02-4030

- 2.5V LED illumination for true tissue color
- Wide angle swivel 2.5X viewing lens permits instrumentation
- 19 viewing lenses from -20 to +20 diopter for better resolution
- "AA" battery handle (batteries not included) with on/off power switch integrated into pocket clip
- 4 year speculum/cannula with internal diameter 2, 3, 4 & 5 mm



Mini Otoscope

02-4010

- 2.5V LED illumination for true tissue color
- Wide angle swivel 2.5X viewing lens permits instrumentation
- "AA" battery handle (batteries not included) with on/off power switch integrated into pocket clip
- 4 year speculum/cannula with internal diameter 2, 3, 4 & 5 mm
- Nylon pouch with belt loop



Mini Ophthalmoscope

02-4020

- 2.5V LED illumination for true tissue color
- 19 viewing lenses from -20 to +20 diopter for better resolution
- "AA" battery handle (batteries not included) with on/off power switch integrated into pocket clip
- Nylon pouch with belt loop



Percussion Hammers

TAYLOR Percussion Hammer

Latex-free

	Weight	Length
02-6010	70g	20 cm/8"



BERLINER Percussion Hammer

Latex-free

	Weight	Length
02-6050	130g	20 cm/8"



DEJERINE Percussion Hammer

Latex-free

	Weight	Length
02-6020	113g	21 cm/8¼"



DEJERINE Percussion Hammer

Latex-free, with needle

	Weight	Length
02-6030	115g	21 cm/8¼"



TROEMNER Percussion Hammer

Latex-free

	Weight	Length
02-6040	210g	22 cm/8¾"



Diagnostic Instruments

Percussion Hammers

BUCK Percussion Hammer

Latex-free, with needle & brush


	Weight	Length
02-6060	110g	18 cm/7"



BABINSKY Percussion Hammer

Latex-free, head Ø 45 mm

	Weight	Length
02-6070	150g	21 cm/8½"



NATIONAL HOSPITAL Percussion Hammer

Latex-free, head Ø 50 mm

	Weight	Length
02-6080	95g	35 cm/14"



Diagnostic Instruments



Crown Click International

